

# Fluorescence Immunoassay Analyzer



Efficient



Convenient



Fast



Accurate

- Dual detection modes
- Multiple parameter detection

.AFI6100B.



## Dual detection modes

## Multiple parameter detection



### Benefits



#### Instant Result

- Deliver results in 3-15 minutes
- Instantly print test report



#### Easy Data Management

- Powerful information management system (direct connection to LIS)



#### Room Temperature Storage

- Test reagents are stable
- Can be stored at room temperature(2-30°C)



#### Reliable Accuracy

- Has Quantitative result
- Good correlation with wet chemistry method results



#### Easy to operation

- 7 inch touch screen
- One-key switch
- Support for standard test & quick test



#### Multiple products offerings including

- Diabetic & Renal Injury, Thyroid Function, Hormone, Cardiac Markers, Inflammation Markers, Tumor Markers, Gastric Function Markers, Rheumatism, Infection Disease, Drug of Abuse, Others

| Test Item | Catalog | Linear range | Specimen | Reading Time | Qualification |
|-----------|---------|--------------|----------|--------------|---------------|
|-----------|---------|--------------|----------|--------------|---------------|

## Tumor Markers

|                                           |             |              |       |         |    |
|-------------------------------------------|-------------|--------------|-------|---------|----|
| <b>NEW</b> Cancer Antigen 125             | CA125-W23F  | 5-500U/mL    | S/P   | 15 min. | /  |
| <b>NEW</b> Cancer Antigen 15-3            | CA15-3-W23F | 1-300U/mL    | S/P   | 15 min. | /  |
| <b>NEW</b> Cancer Antigen 19-9            | CA19-9-W23F | 5-600U/mL    | S/P   | 15 min. | /  |
| Alpha Fetal Protein                       | AFP-W23F    | 2.5-400ng/mL | W/S/P | 15 min. | CE |
| Prostate-specific Antigen                 | PSA-W23F    | 2-100ng/mL   | W/S/P | 15 min. | /  |
| <b>NEW</b> Free Prostate-specific Antigen | FPSA-W23F   | 0.2-30ng/mL  | S/P   | 15 min. | /  |
| Carcino-embryonic Antigen                 | CEA-W23F    | 1-500ng/mL   | W/S/P | 15 min. | CE |
| Ferritin                                  | FRT-W23F    | 10-1000ng/mL | W/S/P | 15 min. | CE |
| Fecal Occult Blood                        | FOB-F23F    | 25-1000ng/mL | F     | 15 min. | CE |

## Diabetic & Renal Injury

|                  |            |             |       |                             |    |
|------------------|------------|-------------|-------|-----------------------------|----|
| HbA1c            | HbA1c-W23F | 4-15%       | W     | Lyses 1 min. & React 5 min. | CE |
| Albumin          | ALB-U23F   | 5-300µg/mL  | Urine | 5 min.                      | CE |
| Cystatin C       | CysC-W23F  | 0.2-10µg/mL | W/S/P | 5 min.                      | CE |
| β2-microglobulin | β2-MG-W23F | 0.3-20µg/mL | W/S/P | 10 min.                     | CE |

## Thyroid Function

|                                    |            |                 |       |                              |    |
|------------------------------------|------------|-----------------|-------|------------------------------|----|
| TT3                                | TT3-W23F   | 0.61-9.22nmol/L | W/S/P | Lyses 1 min. & React 10 min. | CE |
| TT4                                | TT4-W23F   | 10-310nmol/L    | W/S/P | Lyses 1 min. & React 10 min. | CE |
| TSH                                | TSH-W23F   | 0.1-100mIU/L    | W/S/P | 15 min.                      | CE |
| FT3                                | FT3-W23F   | 0.4-50pmol/L    | W/S/P | Lyses 5 min. & React 15 min. | CE |
| FT4                                | FT4-W23F   | 1-100pmol/L     | W/S/P | 15 min.                      | CE |
| <b>NEW</b> Anti-Thyroid peroxidase | A-TPO-W23F | 10-600IU/mL     | W/S/P | 15 min.                      | /  |
| <b>NEW</b> Thyroglobulin           | TG-W23F    | 5-500ng/mL      | W/S/P | 15 min.                      | /  |

## Hormone

|                              |            |                |       |         |    |
|------------------------------|------------|----------------|-------|---------|----|
| β-HCG                        | β-HCG-W23F | 5-200000mIU/mL | W/S/P | 15 min. | CE |
| Anti-mullerian Hormone       | AMH-W23F   | 0.1-16ng/mL    | W/S/P | 15 min. | CE |
| Luteinizing Hormone          | LH-W23F    | 1-100mIU/mL    | W/S/P | 15 min. | CE |
| Follicle-Stimulating Hormone | FSH-W23F   | 1-100mIU/mL    | W/S/P | 15 min. | CE |
| Progesterone                 | Prog-W23F  | 1.4-60ng/mL    | W/S/P | 15 min. | CE |
| Prolactin                    | PRL-W23F   | 1-200ng/mL     | W/S/P | 15 min. | CE |
| Testosterone                 | TES-W23F   | 0.2-15ng/mL    | W/S/P | 15 min. | CE |
| E2                           | E2-W23F    | 9-3000pg/mL    | W/S/P | 15 min. | CE |

## Rheumatism

|                               |          |             |       |         |   |
|-------------------------------|----------|-------------|-------|---------|---|
| <b>NEW</b> Antistreptolysin O | ASO-W23F | 10-800IU/mL | W/S/P | 15 min. | / |
| <b>NEW</b> Rheumatoid factors | RF-W23F  | 5-600IU/mL  | W/S/P | 15 min. | / |

## Inflammation Markers

|                    |          |              |       |         |    |
|--------------------|----------|--------------|-------|---------|----|
| Procalcitonin      | PCT-W23F | 0.1-100ng/mL | W/S/P | 15 min. | CE |
| C-reactive Protein | CRP-W23F | 0.5-200mg/L  | W/S/P | 3 min.  | CE |
| Interleukin-6      | IL6-W23F | 5-3000pg/mL  | W/S/P | 15 min. | CE |
| Serum Amyloid A    | SAA-W23F | 1-300µg/mL   | W/S/P | 15 min. | CE |

| Test Item | Catalog | Linear range | Specimen | Reading Time | Qualification |
|-----------|---------|--------------|----------|--------------|---------------|
|-----------|---------|--------------|----------|--------------|---------------|

## Cardiac Markers

|                                                     |            |                                                            |       |         |    |
|-----------------------------------------------------|------------|------------------------------------------------------------|-------|---------|----|
| NT-ProBNP                                           | NTBNP-W23F | 20-30000pg/mL                                              | W/S/P | 15 min. | CE |
| Cardiac Troponin I                                  | TRO-W23F   | 0.1-40ng/mL                                                | W/S/P | 15 min. | CE |
| Myoglobin                                           | MYO-W23F   | 5-500ng/mL                                                 | W/S/P | 15 min. | CE |
| Creative Kinase MB                                  | CKMB-W23F  | 2-100ng/mL                                                 | W/S/P | 15 min. | CE |
| Cardiac Troponin I/<br>Myoglobin/Creative Kinase MB | CMC-W23F   | cTnI 0.05-40ng/mL<br>Myo 10-800ng/mL<br>CK-MB 0.5-150ng/mL | W/S/P | 15 min. | CE |
| D-Dimer                                             | DDM-W23F   | 0.1-10µg/mL                                                | W/P   | 10 min. | CE |
| cTnT                                                | cTnT-W23F  | 0.03-10ng/mL                                               | W/S/P | 15 min. | CE |

## Gastric Function Markers

|              |            |            |       |         |    |
|--------------|------------|------------|-------|---------|----|
| Propepsin I  | PG I-W23F  | 2-300ng/mL | W/S/P | 15 min. | CE |
| Propepsin II | PG II-W23F | 1-100ng/mL | W/S/P | 15 min. | CE |

## Others

|                        |           |               |       |                              |    |
|------------------------|-----------|---------------|-------|------------------------------|----|
| Vitamin D3             | VD-W23F   | 10-100ng/mL   | W/S/P | Lyses 5 min. & React 10 min. | CE |
| <b>NEW</b> Vitamin B12 | VB12-W23F | 100-2000pg/mL | S/P   | Lyses 13 min & React 10 min. | CE |

## Infection Disease Tests

| Test Item            | Catalog No. | Specimen        | Qualification |
|----------------------|-------------|-----------------|---------------|
| COVID-19 Ag          | COV-S23F    | NP/OP           | CE            |
| COVID-19 Nab         | COV-W23F    | W/S/P           | CE            |
| RSV                  | RSV-S23F    | Nasal Secretion | CE            |
| Strep A              | STA-S23F    | Throat swab     | CE            |
| Astrovirus           | ASV-F23F    | Feces           | CE            |
| Adenovirus           | ADV-F23F    | Feces           | CE            |
| Adenovirus           | ADV-S23F    | Nasal secretion | CE            |
| Influenza A/B        | FLU-S23F    | Nasal secretion | CE            |
| HIV                  | HIV-W23F    | W/S/P           | /             |
| Syphilis             | SYP-W23F    | W/S/P           | CE            |
| HCV                  | HCV-W23F    | W/S/P           | /             |
| HBsAg                | HBsAg-W23F  | W/S/P           | /             |
| C. difficile Toxin A | CDA-F23F    | Feces           | CE            |
| C. difficile Toxin B | CDB-F23F    | Feces           | CE            |
| C. difficile GDH     | GDH-F23F    | Feces           | CE            |
| H.Pylori Ab          | HP-W23F     | W/S/P           | CE            |
| H.Pylori Ag          | HP-F23F     | Feces           | CE            |
| Rotavirus            | ROTA-F23F   | Feces           | CE            |

## Drug of Abuse Tests

| Test Item              | Catalog No. | Specimen | Qualification |
|------------------------|-------------|----------|---------------|
| Amphetamine            | AMP-H23F    |          | CE            |
| Barbiturate            | BAR-H23F    |          | CE            |
| Benzodiazepine         | BZO-H23F    |          | CE            |
| Cocaine                | COC-H23F    |          | CE            |
| Fentanyl               | FYL-H23F    |          | CE            |
| Ketamine               | KET-H23F    |          | CE            |
| K2 Synthetic marijuana | K2-H23F     | Hair     | CE            |
| Marijuana              | THC-H23F    |          | CE            |
| MDMA                   | MDMA-H23F   |          | CE            |
| Methamphetamine        | MET-H23F    |          | CE            |
| Morphine               | MOR-H23F    |          | CE            |
| Tramadol               | TML-H23F    |          | CE            |
| 6-Acetylmorphine       | 6-MAM-H23F  |          | CE            |
| Fentanyl(0.5/1 ng/mL)  | FYL-U23F    | Urine    | /             |



**Assure Tech. (Hangzhou) Co.,Ltd.**

Building 4, No. 1418-50, Moganshan Road,  
Gongshu District, Hangzhou, 310011 Zhejiang, P.R. China



**400-853-5577**